

EAR TOILET PACK






OPERATING MICROSCOPE

Alltton YSX LED Operating Microscope

ALLTTON is proud to offer YSX series of operation microscopes. Precise structure, modern appearance and LED illumination system. Offering excellent parfocal, large field view, precise and clear three dimensional images. The YSX series of operation microscopes have optional working distances and magnification and are suitable for numerous departments including ophthalmology, stomatology, gynaecology, ENT and surgery.

- 1500mm Super long pantograph arm offers a wide operation space
- Four different working distances: 200mm, 250mm, 300mm, 400mm
- Smooth turret 3-step magnification pod
- Field diameter: 7.88mm-81mm
- Magnification: 2.4X-12X
- Green and blue filter selection
- 10W LED illumination
- 20,000 hours LED lifespan
- 50,000 Lux object illumination
- Power supply voltage: 90V/60HZ-250V/50HZ



SUCTION & ACCESSORIES


Tongye Suction Unit

- High vacuum, high flow, oil free vacuum pump
- Antibacterial filter
- Anti-overflow Protection Device
- Lower noise, stepless vacuum adjustment
- Multi-function Deployment Optional
- AC Power Supply



Ear Speculum Gruber

Stainless steel. Ear Gruber speculum are used for opening the ear canal for better visualization during cerumen management.



Suction Controller

A controller is needed to insert the suction tube and controls the level of suction through thumb vent mechanism.




Code	Description	Colour	Length	Material	QTY
NZ3238100	Semi-Disposable Suction Controller, Non Sterile	Black	39mm	ABS	10

Micro Ear Suction Cannula

Suction tubes for suction of fluid and other media from the ear.

- Suction tubes have colour coded hubs for easy identification of tube size
- Cost-efficient and contamination risk free alternative to reusable versions



Code	Description	Outside Diameter	Colour	Length	QTY
NZ5254Bent	Suction Tube 14 Gauge Bent 100 Pack	2.1 mm	Green	100 mm	100/Pack
NZ5250Bent	Suction Tube 16 Gauge Bent 100 Pack	1.6mm	White	100 mm	100/Pack
NZ5251Bent	Suction Tube 18 Gauge Bent 100 Pack	1.2 mm	Pink	100 mm	100/Pack
NZ5252Bent	Suction Tube 20 Gauge Bent 100 Pack	0.9mm	Yellow	100 mm	100/Pack
NZ5253Bent	Suction Tube 22 Gauge Bent 100 Pack	0.7mm	Grey	100 mm	100/Pack
NZ5255Bent	Suction Tube 24 Gauge Bent 100 Pack	0.5mm	Purple	100 mm	100/Pack

MEDICARE REBATE # 41647

FEE \$109.90

Benefit 75%

Benefit 85%

\$82.45

\$93.45

