

# MICRO EAR SUCTION TUBES & CONTROLLERS

## SUCTION CONTROLLER

A controller is needed to insert the suction tube and controls the level of suction through thumb vent mechanism.




Code	Description	Colour	Length	Material	QTY
NZ3238100	Semi-Disposable Suction Controller, Non Sterile	2.1 mm	39mm	ABS	10

## MICRO EAR SUCTION TUBES

Suction tubes for suction of fluid and other media from the ear.

- Suction tubes have colour coded hubs for easy identification of tube size
- Cost-efficient and contamination risk free alternative to reusable versions



Code	Description	Outside Diameter	Colour	Length	QTY
NZ5254Bent	Suction Tube 14 Gauge Bent 100 Pack	2.1 mm	Green	100 mm	100/Pack
NZ5250Bent	Suction Tube 16 Gauge Bent 100 Pack	1.6mm	White	100 mm	100/Pack
NZ5251Bent	Suction Tube 18 Gauge Bent 100 Pack	1.2 mm	Pink	100 mm	100/Pack
NZ5252Bent	Suction Tube 20 Gauge Bent 100 Pack	0.9mm	Yellow	100 mm	100/Pack
NZ5253Bent	Suction Tube 22 Gauge Bent 100 Pack	0.7mm	Grey	100 mm	100/Pack
NZ5255Bent	Suction Tube 24 Gauge Bent 100 Pack	0.5mm	Purple	100 mm	100/Pack