

ssaustralia HEALTHCARE
SUPPLIES

Does More For you

Furniture & Equipment

PRODUCT GUIDE

1800 777 518 | WWW.SSSAUSTRALIA.COM.AU

SSS AUSTRALIA | Supporting Quality Care

SSS Australia is committed to supporting you in providing the highest quality of care for your patients. Our extensive range of furniture and equipment ensures we have the product solution for your facility. We choose to supply high quality, innovative products to help you provide the highest quality care for your patients.

For more information on any of the products in this booklet please do not hesitate to contact your personal account manager or our expert customer support team.

National Representatives

Contact our friendly customer support team to request a quote

Customer Support

1800 777 518

sales@sssaustralia.com.au

Aged Care Representatives

Contact your Aged Care personal account manager to request a quote:

Shannon Jennings

Brisbane South, Ipswich, Toowoomba & South Australia

0419 172 150

shannon.jennings@sssaustralia.com.au

General Healthcare Representatives

Contact your General Healthcare personal account manager to request a quote:

Chris Mitchell

Brisbane & Southside, Ipswich to Gatton, Mackay to Ayr

0448 031 326

chris.mitchell@sssaustralia.com.au

Contact your Regional business development manager to request a quote:

Amaan Vaddadi

Victoria

0400 121 010

amaan.vaddadi@sssaustralia.com.au

Contact your Aged Care personal account manager to request a quote:

Fiona Carey

Brisbane North, Wide Bay, Capricorn & Western Australia

0427 361 761

fiona.carey@sssaustralia.com.au

Contact your General Healthcare personal account manager to request a quote:

Michelle Warnecke

Brisbane North, Brisbane CBD and Far North

0419 289 852

michelle.warnecke@sssaustralia.com.au

Contact your Regional business development manager to request a quote:

Jeff Osbourne

Greater Sydney & New South Wales

0428 759 080

jeff.osbourne@sssaustralia.com.au

Contact your Aged Care personal account manager to request a quote:

Paul Cooper

Sunshine Coast, North QLD & Far North QLD

0482 629 882

paul.cooper@sssaustralia.com.au

Contact your General Healthcare personal account manager to request a quote:

Chris Wands

Sunshine Coast, Wide Bay & Capricorn

0408 061 221

chris.wands@sssaustralia.com.au

Contact your Hospital personal account manager to request a quote:

Doug McDonnell

National

0473 871 135

doug.mcdonnell@sssaustralia.com.au

Contact your Aged Care personal account manager to request a quote:

Carol Edmond

Logan, Gold Coast, Northern Rivers, West Brisbane & Northern Territory

0408 181 938

carol.edmond@sssaustralia.com.au

Contact your General Healthcare personal account manager to request a quote:

Melissa Kneale

Gold Coast, Toowoomba & Northern Rivers

0407 031 599

nick.iles@sssaustralia.com.au

Does More For you

TABLE OF CONTENTS

DAILY LIVING AIDS

4

- Bathroom Aids
- Mobility Aids
- Warmth & Comfort
- Mealtime Aids
- T.E.N.S Therapy
- Nursing Aids
- Crash Mats

MANUAL HANDLING & SAFETY

11

- Slide Sheets
- Patient Slide
- Sit & Slide Pad
- Resident Transfer Belt
- Emergency Resident Transfer
- Hoist Slings
- Hoists
- Patient Transfer Aids
- Safety Socks
- Nurse Call Bed Alarms

MEDICATION & VACCINE FRIDGES

16

- Vaccine Fridges

FURNITURE & TROLLEYS

17

- Beds
- Bedroom Accessories
- Trolleys
- Stainless Steel Trolleys
- Oxygen Trolleys
- Stools
- Arm Chairs

PRESSURE MANAGEMENT

22

- Mattresses
- Mattresses with Pressure Reducing Insert
- Alternating Mattresses
- Mattresses Underlay
- Pressure Cushions
- Body Protection

LINEN & LAUNDRY

32

- Sheets & Pillowcases
- Blankets
- Bath Towels & Face Washers
- Pillows
- Linen Protection
- Laundry Bags
- Linen Skips

MEDICATION ADMINISTRATION

36

- Syringe Pumps
- Medication Rounds

DIAGNOSTIC MEASURE

37

- Blood Pressure - Automatic
- Blood Pressure - Manual
- Vital Signs Monitoring
- Otoscopes
- Thermometers
- Pulse Oximeters
- Bladder Scanners
- Stethoscopes
- Scales
- Exam Lights
- Point of Care Testing
- Anaemia Screening
- Glucometers
- Stand Alone ECG
- Oxygen Concentrators

RESUSCITATION & CRITICAL CARE

46

- Suction Pumps
- Nebuliser Therapy
- Defibrillators

BATHROOM AIDS



Aspire Shower Commodes

Easy-slide weight bearing platform footrest slides completely away for transfers.

- 4x total locking stainless steel castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 8° backrest angle ensures patient comfort
- Improved toilet clearance

Item #	Max Weight	Seat Height	Seat Width	UOM
1376160	200kg	590mm	460mm	Each
1375989	200kg	590mm	530mm	Each
1335162	400kg	590mm	600mm	Each



Aspire Tilt-in-Space Shower Commode

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Height adjustable, swing away footrests with calf pads for optimal positioning, support and ease of transfers
- Tool-free height adjustment to 4 settings allow the commode to be easily adapted to varying needs

Item #	Max Weight	Seat Height	Seat Width	UOM
1368540	250kg	460-630mm	460mm	Each
1408356			530mm	Each



Aspire Padded Commode Seats

- Moulded and seamless soft PU seats designed for comfort and durability
- Choose the size appropriately to fit your Aspire Shower Commode
- Easy clean, fully sealed surface for infection control
- Adaptable to most shower commode brands

Item #	Spec	Fits	UOM
1376268	Grey Open Front 460mm	1376160, 1368540	Each
1404857	Grey Closed Front 460mm	1376160, 1368540	Each
1375997	Grey Open Front 530mm	1375989	Each
1335170	Grey Open Front 600mm	1335162	Each
1396250	Bowl and Pan Carrier	1376160, 1368540, 1375989, 1335162	Each

Smik Shower Commode Sliding Footrest & Open Seat

- Can be positioned over the toilet for toileting
- Swingback arms assist side transfers
- Pan carrier included
- Padded seamless commode seat included
- Patient locking arms
- Weight bearing sliding footrest
- Slide in foot plate
- PVC slatted backrest

Item #	Seat Width	UOM
1454854	550mm	Each
1386867	450mm	Each





Alphacare Lopital Reflex Electric Shower Commode

- Foldable swivel armrests with front over safety locking arms
- Electrically operated height-adjustable settings (variable)
- Easy to recharge with battery charger provided
- Height-adjustable legrests and calf supports
- Ergonomically designed backrest and seat
- Carrier for stainless steel commode pan
- Suitable for use in the shower and toilet
- Easy-to-use touch controls on backrest
- Electrically operated tilt mechanism
- Also available with Hand-control
- Height-adjustable headrest
- 4x braking castors



BULK SHIPMENT

Item #	Max Weight	Seat Height	Seat Width	UOM
1435310	135kg	490 - 1001mm	480mm	Each

Aquarius Mobile Shower Trolley



- Quality 125mm twin wheel castors with total lock and steer
- Liner inner (min) base size - 560 x 1720 x 180mm deep
- 795 W x 1950mm Long overall size including bumpers
- Powder coated stainless steel dropsides and endrails
- Adjustable end rails with 120mm extra length
- 135mm lowest clearance under the base
- Economical hydraulic trolley
- Made in Australia
- Hose Length 1100mm



BULK SHIPMENT

Item #	Max Weight	Height Range	Overall Size	UOM
1373250	180kg	680-880mm	795x1950mm	Each



Etac Swift Mobil-2 Tilt-in-Space Shower Commode

Etac Swift Mobil Tilt-2 with XL Back is designed with both the carer and user's needs in mind. The wider backrest allows for an additional 6cm between armrests to accommodate larger users. It handles the multitude of challenges faced by caregivers in their daily work and provides the comfort and safety users need to feel secure.

Item #	Max Weight	Height Range	UOM
1481843	160kg	680-880mm	Each



Etac Clean Mobile Shower Commode 49cm

The Etac Clean Mobile Shower Chair is an award-winning commode designed to be secure and comfortable for the user. To increase durability the Shower Chair is constructed without holes or welding seams to reduce corrosion. This Shower Chair has a functional push handle which makes it very easy to manoeuvre and four locking castors to keep the user safe when stationary.

Item #	1481851	Each



BetterLiving Mobile Shower Commode

The BetterLiving Mobile Commode is designed for aged care, community and home care environments. The frame and its major components are constructed from corrosion resistant aluminium and powder coated for extra protection. All accessories, including pins, nuts and bolts are constructed with corrosion-resistant materials to help ensure the strength and longevity of the product when used as intended.

Item #	1481878	Each

BATHROOM AIDS



Shower Chairs

- Comfortable and hygienic plastic moulded seat
- Gentle slope allows water to easily flow away



Item #	Max Weight	Seat Height	Seat Width	UOM
1292633	110kg	390-540mm	440mm	Each
1064094	125kg	420-550mm	440mm	Each
1292641	125kg	420-570mm	440mm	Each
1426508	200kg	375-550mm	440mm	Each



Over the Toilet Shower Chairs

This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non-slip rubber feet. Comes complete with a seat and splashguard.



Item #	Max Weight	Seat Height	Seat Width	UOM
1175280	125kg	440-490mm	370mm	Each
1214807	135kg	470-590mm	500mm	Each
1214790	210kg	515mm	650mm	Each



Bariatric All-In-One Over Toilet & Shower Chair Commode

- Comes with square bowl with easy front access to remove bowl
- The seat lifts up and comes with a seat cover
- Removable backrest adds stability when used as a bedside commode

Item #	Max Weight	Seat Height	Seat Depth	UOM
1286250	295kg	465-570mm	450mm	Each



Aspire Deluxe Bedside Commode

- Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle



Item #	Colour	UOM	Max Weight	Seat Height	Seat Depth	Seat Width
1408364	Slate Blue	Each	160kg	410-565mm	450mm	460mm
1434810	Champagne	Each	160kg	410-565mm	450mm	460mm

BATHROOM AIDS



Toilet Seat Reducer

- Designed for thin frail and bony patients who may slip down into the toilet bowl
- Use the Seat Reducer on commode chairs, shower chairs or straight on a toilet seat
- Soft and comfortable foam that will reduce the gap of the hole in the toilet seat

Item #	Hole Size	UOM
1289616	205x80mm	Each



Pelican Shower Chair Cape

- Gives dignity when taking patients to and from the shower room
- Easy to put on by simply putting it over the patient's head
- Polycotton material that will dry quick

Item #	1028169	Each



Betterliving Cognitive Assistance Toilet Seat Blue

The BetterLiving Cognitive Assistance Toilet Seat is vibrant blue in colour to provide assistance for those with vision impairment or cognitive difficulties.

The toilet seat is made of quality Duroplast material and features stainless steel hardware which can be attached to the toilet pan via top or bottom fixing method. Its universal design means it is compatible with most toilets.

Easy to install, the toilet seat is designed with a maximum user weight of 175 kilograms.

Item #	Max Weight	UOM
1492691	175kg	Each

MOBILITY AIDS



Gutter Frame Rollator Walker

The Days Gutter Walker with Handbrakes has a new improved design with thicker tubing and larger wheels. Designed to keep the user in an upright position for gait training. Folds for storage. Forearm and hand supports are fully adjustable.

- Between Wheels: 580mm

Item #	Max Weight	Height	Width	UOM
1304534	150kg	630mm	630mm	Each



Seat Walkers with Basket & Handbrakes

Item #	Colour	Max Weight	Height	UOM
1250017	Black	150kg	830-930mm	Each
1396781	Red	150kg	830-930mm	Each
1396773	Blue	150kg	830-930mm	Each
1498592	Beige			Each
1503855	Black	150kg		Each



MOBILITY AIDS



Days Link Wheelchair

- Lightweight, adaptable
- Self Propelled
- Practical

Item #	Spec	UOM
1292986	20" 150kg	Each
1284290	18" 130kg	Each



BULK SHIPMENT



Days Whirl Wheelchair

- 18" wide seat with 120kg user capacity
- Swing away and removable leg rests
- Puncture proof solid tyres
- Removable armrests
- Robust steel frame

Item #	Spec	UOM
1279130	Self Propelled	Each
1279149	Transit Attendant Propelled	Each



BULK SHIPMENT



Days Swift Wheelchair

- Fixed Arms Swing away and removable leg rests
- 18" wide seat with 115kg user capacity
- Puncture proof solid tyres
- Good manoeuvrability
- Robust steel frame
- Ultralight 7.8kg

Item #	Spec	UOM
1324810	Self Propelled	Each
1331957	Transit Attendant Propelled	Each



BULK SHIPMENT



Sentra Bariatric Wheelchair

- Swing back removable armrests
- Sturdy, carbon steel frame
- Simple folding mechanism
- 200kg max user weight
- Swing in/out leg rests

Item #	Spec	UOM
1274420	24" 600mm	Each
1273183	22" 550mm	Each



BULK SHIPMENT



ID Soft Tilt-In-Space Reclining Wheelchair

- Sturdy steel frame providing extra support and stability in any position for added security
- Infinite position adjustment of the backrest and seat
- Larger wheels for a smoother ride
- 130kg Maximum user weight

Item #	Spec	UOM
1297562	400mm 18"	Each
1027635	Footplate Covers	Pair



BULK SHIPMENT

WARMTH & COMFORT



Beurer Electric Heat Warmers

- Cosy and soft – cover made from microfibre
- Overheating protection
- Automatic-off function after approx. 90 minutes
- 3 temperature settings
- Machine-washable at 30°C

Item #	Spec	UOM
1414909	44x33cm	Each



Sheepskin Underblanket

- Natural fibres for better temperature balance
- Withstands high washing temperatures
- Urine resistant

Item #	Spec	UOM
1172672	30mm	Each

MEALTIME AIDS



Rolyan Manoy Built Up Plates

- Plates with sloping bottom are ideal for people who have the use of only one hand
- Dishwasher safe
- Not microwave safe

Item #	Size	UOM
1053368	Small 220x180mm	Each
1053359	Large 280x200mm	Each
1028347	Plate Surround	Each



Felt Clothing Protectors

- Highly absorbent with impermeable barrier backing
- Dignified protection for personal clothing
- Adjustable snap closure

Item #	Colour	Spec	UOM
1221490	Azure Blue	450x450mm	Each
1289659	Blue Plaid	900x450mm	Each
1287028	Maroon Plaid	900x450mm	Each



White Terry Towelling Clothing Protector

Item #	Spec	UOM
1215368	Pop Over Style with Yellow Binding	Each



Rolyan Supagrip Cutlery Range

Item #	Spec	UOM
1184889	Knife	Each
1184898	Fork	Each
1184904	Tablespoon	Each
1395319	Teaspoon	Each

T.E.N.S THERAPY



Omron® Deluxe HV-F128

The Omron HV-F128 is the deluxe T.E.N.S device from Omron. With twelve automatic therapy modes the Omron HV-F128 can easily be set to offer either pain relief or muscle stimulation at the touch of a button.

- Provides relief of pain
- Suitable for management of muscular injuries
- High frequency stimulation to relieve pain quickly
- Single channel device

Item #	1242006	UOM	Each
--------	---------	-----	------



Omron® Standard HV-F127

Compact in its design the Omron HV-F127 is a simple to use T.E.N.S device. With one touch operation you can easily select one of the nine automatic therapy modes.

- High frequency stimulation to relieve pain quickly
- Suitable for management of muscular injuries
- Provides relief of pain
- Single channel device

Item #	1241591	UOM	Each
--------	---------	-----	------



ProTens Dual Channel T.E.N.S Machine Kit

Kit contains:

- One dual-channel T.E.N.S stimulator
- Carrying case
- Lead wires
- Operating manual
- 9V batteries,
- Pack of four electrodes 50mm

Item #	1280813	UOM	Each
--------	---------	-----	------



Omron® T.E.N.S Electrodes

- Suitable for use with Omron HV-F128 & HVF127
- No messy gels required
- Long Life Pads

Item #	1170373	UOM	Pack 2
--------	---------	-----	--------



Pals T.E.N.S Electrodes

To suit Protensens machine

Item #	Spec	UOM
1326163	Valutrode 50x50mm Square	Pack 4
1252316	Platinum 50x50mm Square	Pack 4
1054365	Platinum 50mm Round	Pack 4
1064058	Platinum 50x90mm Rectangle	Pack 4

NURSING AIDS



Pelican™ Smoking Aprons

- Designed to protect the smoker from burns if they drop a cigarette
- Made from kevlar and nomex (used by firefighters and racing car drivers)
- Can be washed without losing its fire retarding properties

Item #	Spec	UOM
1258919	Without Side Flaps	Each
1287354	With Side Flaps	Each

CRASH MATS



PremiumLift Crash Mat

- Helps reduce impact forces by up to 90%
- Provides full coverage of the bed exit area
- 3cm thick foam creates low profile edge
- Tri-fold design for easy storage
- TGA Approved
- Wipe clean vinyl cover

Item #	Spec	UOM
1453726	1780x965x30mm	Each



HealthSaver Bedside Crash Mat

- Made from high-impact foam with an easy-to-clean vinyl cover
- Non-slip material ensures that the floor mat remains firmly in place
- Sloping edges make it easy to manoeuvre a wheelchair over the fall mat

Item #	Spec	UOM
1325485	1820x1000x50mm	Each

SLIDE SHEETS



Silicone Turning Slide Sheets

- These special slippery sheets are used for turning residents in bed and providing easier manual handling
- Made from washable, drip-dry, light-weight fabric that is silicone coated

Item #	Spec	Size	Colour	UOM
1240413	Roll	50 x 1.5m	Blue	Roll
1274037	Double	2 x 1.5m	Red	Each
1396307	Double	2 x 1.5m	Green	Each



Slide Sheet Holder

Item #	Colour	UOM
1362720	Blue	Each
1334638	Red	Each
1293311	Colour Varies	Each



WendyLett Slide Sheets

The WendyLett base sheet is a unique sliding sheet that appears to be a white and grey sheet with narrow white stripes, and is available in many different sizes. The sliding takes place sideways along the direction of the stripes.

Item #	Spec	Size	UOM
1331353	4 Way Draw King Single	2000x2000mm	Each
1393650	4 Way Draw Single	1400x2000mm	Each
1331396	Base Sheet King Single	1050x2000mm	Each
1331361	Base Sheet Single	900x2000mm	Each

PATIENT SLIDE



Patslide® All Purpose Resident Slide

- An efficient, cost effective and comfortable method for general resident transfers
- Light weight, radiolucent, easy-to-clean and maintenance free
- Safe and easy to use
- Static resistive

Item #	Spec	UOM
1078788	1525x635mm	Each

SIT & SLIDE PAD



Pelican™ Sit & Slide Chair Pad

- Allows resident to sit down at the front of a chair and then slide hips back into the chair
- One or two nurses can give assistance by pulling the pad back with the lifting loops
- The resident can also be helped to stand by carer pulling the webbing handles

Item #	Spec	UOM
1125841		Each

RESIDENT TRANSFER BELT



Pelican™ Soft Transfer Belt

- Four large padded vertical handles
- Easily adjusted for size to suit a variety of residents
- The soft padded belt is slip resistant

Item #	Spec	UOM
1164358	800x1600 mm	Each

EMERGENCY RESIDENT TRANSFER



AlbacMat Rescue Mat

- Designed with the safety of the rescuer and the patient in mind
- Able to travel easily over most surfaces - including stairs, carpet, and concrete - allows for improved efficiency when evacuating a facility
- Evacuations can be performed by one person per patient in most situations which is essential when time is critical
- Its compact, self-storage, and lightweight design means it can be readily available in varying locations for instant usage

Item #	Spec	Load	UOM
1471805	Standard 600x1850mm	<100kg	Each
1471813	Large 900x1850mm	>100kg	Each



AlbacMat Holder

Holds four mats vertically and allows for quick access in an emergency.

Item #	1471821	Each



ResQmat Pelican Evacuation Mat Padded

- Completely encloses the legs and arms. This prevents people from reaching out whilst being rescued, and possibly impeding the rescue
- Full length mini mattress protects the body when going down stairs
- Padded head pillow protects the head when going down stairs
- Strong, thick Hook & Loop Tape straps
- Automatically fits to size

Item #	Max Weight	UOM
1302310	300kg	Each



Pelican™ Patient Rescue Sheet

- The mattress and the resident can then be lowered onto the floor for very fast and safe evacuation, without staff having to carry them
- In case of fire or emergency, the webbing straps are pulled out and quickly done up around the resident
- Kept in position under the mattress with strong elastic corner straps
- Daily bed making is as normal

Item #	Spec	Max Weight	UOM
1027949	Rescue Sheet for Vinyl Only	300kg	Each
1249740	Rescue Sheet for Vinyl & Carpet	300kg	Each




HOIST SLINGS



Pelican™ Premier Lifting Sling

- Full body sling with soft padded split legs and padded head support
- Designed to give a variety of lifting positions to enable residents to be lifted in a reclined mode or in a sitting upright position

Item #	Spec	UOM
1027840	White Nylon	Each
1143385	Blue Mesh	Each



AidaCare Aspire Deluxe Swings

- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting

Item #	Spec	Max Load	UOM
1415409	General Purpose Sling with Head Support	200kg	Each
1425587	General Purpose Sling with Head Support	200kg	Each
1425609	General Purpose Sling with Head Support	300kg	Each
1425595	General Purpose Sling with Head Support	300kg	Each
1415220	Hygiene Access Sling with Head Support	200kg	Each
1415204	Hygiene Access Sling with Head Support	200kg	Each
1415212	Hygiene Access Sling with Head Support	300kg	Each
1415239	Hygiene Access Sling with Head Support	300kg	Each
1425560	Padded Standing Sling	200kg	Each
1425552	Padded Standing Sling	300kg	Each

HOISTS



Aspire 200 Steel Lifter

- Designed to meet the rigorous demands of professional and facility environments
- Ergonomic handles provide multiple grip points to suit carers of different heights and strengths
- Powered leg function opens and closes legs for lifting stability or accessibility
- Strong powder coated steel frame with industry leading Linak electronics
- Critical safety features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- Warranty - Frame 5 Years, Electronics 3 Years

Item #	Leg Gait	Lift Range	Max Load	UOM
1373056	580-900mm	535-1725mm	200kg	Each



Aspire 200 Steel Standing Lifter

- Bullhorn Lifting Handles for patient stability and 6 sling attachment points
- Ergonomic Handles Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function Opens and closes the legs for lifting stability or accessibility
- Warranty - Frame 5 Years, Electronics 3 Years, Battery & Handset 1 Year

Item #	Leg Gait	Lift Range	Max Load	UOM
1391054	575-936mm	850-1740mm	200kg	Each



Aspire A205 Aluminium Patient Lifter

- The Aidacare Aspire A205 is a compact and versatile aluminium mobile lifter with powered legs that is ideal for home, community and facility environments
- Transportation or Storage: Is quick and convenient with tool free disassembly
- Ergonomic Handset: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Adjustable Height: Three adjustable height settings provides the caregiver with the flexibility to adjust the height range to suit the specific lift and patient scenario

Item #	Leg Gait	Lift Range	Max Load	UOM
1422766	547-1015mm	395-1805mm	205kg	Each



HOISTS



Aspire Inline Weigh Scale

- Designed to be fitted on Aidacare Aspire Patient Lifters
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display

Item #	Max Weight	UOM
1415727	320kg	Each



Etac Turner PRO Standing Aid

The Etac Turner PRO Patient Transfer Aid provides support for a patient during transfers between seated positions and provides ergonomically correct positioning for the caregiver. Using this product reduces the risk of falls and enables better communication between the patient and caregiver.

The orange handle provides an intense, high contrasting colour that can help users with impaired vision or dementia.

Etac Turner Pro is compact and easy to grip, and due to its low weight of only 7.5 kg, it is easy to transport.

Item #	Max Weight	UOM
1481886	200kg	Each

PATIENT TRANSFER AIDS



Aspire Star Patient Mover

The Aspire Star Patient Mover is an ideal alternative to a patient lifter for users with appropriate sitting balance, who require transfers over greater distances.

- Height adjustable knee support ensures correct positioning for a larger patient height range
- Low footplate height makes patient entry and exit easier
- Wide seat entry suits a larger range of body shapes

Item #	Total Width	Total Height	Max Load	UOM
1461575	615-900mm	1035mm	185kg	Each

SAFETY SOCKS



SafeStep Ankle Height Safety Socks

- Non-slip grip spots on top and bottom of socks
- Bright orange for easy identification of residents at risk


Item #	Size	Colour	UOM
1320157	Regular Sml-Med (Sizes 2-6)	Orange	Pair
1320173	Large (Sizes 7-11)	Orange	Pair
1320181	X-Large (Sizes 12-15)	Orange	Pair
1459422	X-Large Wide (Size 12-15)	Orange	P12
1323059	Small (Sml-Med Sizes 5-7)	Black	Pair
1323040	Medium (Med-Lge Sizes 7-11)	Black	Pair
1472356	X-Large Wide (Size 12-15)	Red	C480
1472364	Bariatric	Red	C480
1459414	Bariatric	Orange	P12



Foot Slipper Safety Socks

- Double sided grip
- Long

Item #	Size	Colour	UOM
1331191	Lge-XL	Red	Pair
1331213	Lge-XL	Dark Green	Pair



Out of bed alarms are an inexpensive and highly effective method of reducing falls by monitoring the out of bed movements of at risk residents.

The Out of Bed Sensor Mat can be easily fitted to your existing nurse call system. The Sensor Mat interface control replaces your existing nurse call pendant (nurse call push button on a cord) with a new pendant that also connects the Out of Bed sensor mat to your nurse call system. The new pendant & interface control feature an on/off switch for your convenience.

Sensor Beam



- The ultimate early warning for all preventions.
- Utilises motion detection
- Creates an invisible beam above or beside a bed/chair or across a doorway
- Generates a nurse call alarm when beam is interrupted
- Connects to existing nurse call system
- Great for residents who have become aware of floor mats
- Will not false alarm
- Transferable from room to room

Item #	Spec	UOM
1281321	Nurse Call Alarm Sensor Beam w 3 Metre Call Bell Pole & Base	Each
1369121	Nurse Call Alarm Sensor Beam Call Bell with Wall Mount Bracket	Each
1373153	Nurse Call Alarm Sensor Beam with Small Pin RCA	Each
1402005	Sensor Beam w/ Phone Style Call Plug, Bell, Pendant, Pole & Base	Each
1527762	Smart Caregiver Cordless Floor Mat 6096mm x 9144mm	Each
1527770	Smart Caregiver Cordless Floor Mat 6096mm x 1219mm	

Sensor Mats



Item #	Spec	UOM
1415506	Stand Alone Out of Bed Sensor Mat Grey 900 x 600m	Each
1335413	Stand Alone Out of Bed Sensor Mat Grey 1200 x 600mm	Each
Both 1415506 and 1335413 require Interface Control		

Interface Controls



Item #	Spec	UOM
1281305	Interface Control & Pendant A Stereo (2 Rings)	Each
1332678	Interface Control & Pendant A Mono (1 Ring) No Nurse Call	Each
1281410	Interface Control & Pendant B 1 Pin	Each

Please speak to your personal account manager about connection requirements for your facility.

VACCINE FRIDGES

AQ MediSafe Plus Vaccine Fridges

Medisafe Plus Fridges are some of the most feature-packed, high quality models of vaccine storage units available in Australia. Designed and manufactured in Australia, they are purpose built for vaccine storage in line with the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines. Medisafe Plus fridges are QCPP accredited and TGA approved, giving you confidence that your vaccines will be stored securely.

All Medisafe Plus units are uniquely equipped with two temperature probes for rigorous monitoring of both cabinet temperature and a simulated vaccine temperature (when using the AQ Box cloud temperature monitoring system). This yields the most accurate measurements possible and is a unique feature in the market of vaccine protection.

The AQ Box comes standard with all Medisafe Plus fridges, and includes 12 months' free access to the portal. The AQ Box acts as an automated data logger which allows you to view and print temperature charts with Minimum and Maximum temperatures, as well as real time data from any web-enabled smartphone, tablet, or computer. If there is a temperature breach, loss of connectivity or power outage, alerts are sent via SMS and email. Medisafe Plus fridges also come standard with a BMS outlet for hospital remote monitoring.



- Digital temperature control system preconfigured to hold temperature from +2°C to +8°C in line with national 'Strive for 5' storage guidelines
- AQ Box cloud temperature monitoring with SMS and email alerts of temperature breach, loss of connectivity or power outage
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature
- QCPP cold chain accredited and TGA approved
- Audible and visual alarms
- Adjustable shelves
- Lockable door
- Low energy consumption
- Low noise output
- Dedicated technical service line
- 24 hour product support
- 2 year warranty

Accessories

Item #	Capacity	UOM
1279211	80L	Each
1244486	306L	Each
1244495	381L	Each
1415808	1000L	Each



Item #	Spec	UOM
1245400	Spare Shelf for 306L & 381L Vaccine Fridge	Each
1244501	Auto Dialler Fridge Option	Each
1385917	AQ Box Cloud Based Temp Monitoring Device	Each
1329839	Logtag Software & Docking Station Only	Each
1329847	Logtag Temperature Data Logger with Probe	Each



AQ MediSafe Fridge UPS & Battery Pack

The UPS and battery pack provides up to 8 hours run time when attached to a medium sized refrigerated cabinet (ave 300L). Units can be combined with an additional battery pack unit (2+) to increase runtimes of up to 7-8 additional hours per unit. *Run times are based on a new unit and will reduce with usage.

- 3 year warranty

Item #	1436813	Each
--------	---------	------






BEDS

CodaCare

CodaCare is 100% Australian owned. Their story starts way back in 1967, and to this day they continue to build on their heritage of user-centric design, quality local manufacturing, importing and outstanding client service. As leading hospital & aged care residential furniture suppliers, they believe in providing long term value for money for Australia's hospitals, rehabilitation centres, aged care facilities and other healthcare organisations.

They use a range of materials to produce furniture which includes hospital grade compact laminate, a marine grade board that comes with a 10 year warranty. HPL is impact and scratch resistant and does not fade over time.

When you buy Australian made furniture, not only are you supporting local manufacturers and retailers, you are choosing furniture made from quality and sustainable materials which are manufactured under the guise of quality management system such as ISO 9001 and Australian health and safety requirements.



Premium Lift Ultra Low Bed Beech

The PremiumLift Ultra Floorline Bed is the ultimate in low level beds, designed in Germany to the higher quality specification. It offers a variety of functions including Trendelenburg and cardiac chair position allowing for excellent care to be provided. It is an optimal floor level bed solution for use both for home care and aged care facilities. Its minimum height from the floor is just 6.7cm ensuring complete patient safety. Should the user be inclined to roll whilst in bed, and potentially fall, the bed reduces impact and allows the user to safely roll onto the floor, or a fall mat.

Item #	Size	Height	Max Weight	UOM
1427970	Single	70-640mm	185kg	Each
1427962	King Single	70-640mm	185kg	Each



Premium Lift Multi-Flex Mattress

- Castellated gel-infused top layer
- High-Density foam base
- Hinged firm strengthened transfer sides
- 4-Way stretch antimicrobial cover
- Waterproof and resistant to patient body fluids
- Resistance to cleaning agents
- Stain resistant
- Heavy-Duty PVC base
- 360° / Two-Sided integrated zipper
- High-frequency sealed & welded seams
- Waterfall fluid penetration flap

Item #	Size	Spec	Max Weight	UOM
1508814	Single	1980x880x150mm	180kg	Each
1508806	King Single	1980x1060x150mm	140kg	Each

BEDROOM ACCESSORIES



ThermoForm Over Bed Table

This overbed table features a solid "C" base frame and flat top, ideal for use in the hospital and aged care environments. Supplied with a standard willuna white (beige) laminate top with a vinyl wrap around the outside, or with the deluxe "thermoform" willuna white (beige) top which is heat resistant and ideal for hospitals higher requirement needs.

- Thermo heat resistant top with raised edge
- Under-table strata for improved strength
- Lockable castors

Item #	Colour	UOM
1283324	Willuna White	Each



Mobile Stainless Steel I.V. Pole

- Five star solid steel base
- Finished in hammertone gunmetal grey powder coat with black 50mm twin wheel castors
- Chrome plated stem with key wheel for height and tension adjustment

Item #	Colour	UOM
1087357		Each



CodaCare Bedside Locker 3 Draw

- Lockable top drawer providing secure storage for medication and/or valuable personal belongings
- Large, easy-to-grab handles
- Designed with the aged care environment in mind – strong and durable
- Fitted with 4 castors

Item #	Size	UOM
1264952	800x450x500mm	Each



Bed Cradle Blanket Support Epoxy

- Designed to fit under the mattress to keep linen clear of the feet and legs

Item #	Colour	UOM
1096042		Each



Bed Rail Protectors

Item #	Size	Spec	UOM
1373617	500x1660mm	Slip Over	Pair
1027377	500x1660mm	Velcro on Waterproof	Pair

TROLLEYS



Mobile Medicine Trolley 5 Drawer w/ Accessories

Features:

- 304 Stainless Steel Wire Basket
- Luxurious push handle
- Central or Seals Lock
- Accessories included
- Soft closing drawers
- Castors with brakes
- 2 year Warranty
- Extension table
- Non slip mat
- Waste bin

Drawer Sizes	4 drawers x 1250mm	1 drawer x 260mm
--------------	--------------------	------------------

Item #	Colour	Dimensions	Weight	UOM
1337521	White	670x500x1080mm	42kg	Each
1448951	Side Basket To Suit Med Trolley			





Viva Narrow Cart Aluminium 34"

The Aluminium Narrow Cart has many of the same features as the traditional Viva Narrow cart, yet with sleek lines, contemporary handles and a lightweight body for a greater ease of use.

- 4 x 100mm twin swivel castors; two locking
- ABS top with slide out work surface
- Central key lock control
- Push handle
- Soft close drawers

Drawer Sizes	2 drawers x 760mm	2 drawers x 1560mm	1 drawer x 2340mm
--------------	-------------------	--------------------	-------------------

Item #	Colour	Dimensions	Weight	UOM
1451154	Lavender	510x5420x9770mm	49kg	Each
1453076	White	510x5420x9770mm	49kg	Each
1450735	Light Blue	510x5420x9770mm	49kg	Each



Viva Isolation Cart Yellow

The bright yellow colour of the Viva Isolation Cart is synonymous with the high risk infection control areas in the facility. Safe and dust free storage of clean gowns, caps, gloves etc provide peace of mind for staff and visitors.

- 4 x 125mm Twin swivel castors; two locking, one directional
- ABS top with slide-out work surface
- Central key lock control
- Divider in two drawers
- Soft close drawers
- Push handle

Isolation accessory package includes:

- Waste container set and lid
- Multi-storage basket
- Glove Dispenser

Item #	Spec	Dimensions	UOM
1453254	GC1040 3 Drawer	690x502x930mm	Each
1452959	GC1050 3 Drawer	690x520x1010mm	Each
1453203	GC1060 4 Drawer	690x520x10850mm	Each
1452967	Accessory Pack		Each



Emergency Trolley 5 Drawer

This rugged portable crash cart is constructed with stainless steel & an ABS modular surface. With a heavy load rating and complete with accessories it's a must for many healthcare facilities.

- Noise resistant castors with brakes
- 304 stainless steel wire basket
- Luxurious push handle
- External power holder
- Stainless steel IV pole
- ABS drawer dividers
- Central or seals lock
- Oxygen tank holder
- Soft closing drawer
- Defibrillator Shelf
- Non slip mat
- CPR board
- Waste bin

Drawer Sizes	3 drawers x 110mm	1 drawer x 1350mm	1 drawer x 1950mm
--------------	-------------------	-------------------	-------------------

Item #	Colour	Dimensions	Weight	UOM
1324926	Red	650x480x980mm	55kg	Each



STAINLESS STEEL TROLLEYS

Econo Trolleys

- Welded stainless steel construction with safety fold to edges
- Smooth 95mm noise-free castors with front brakes
- Magnetised drawers to prevent opening
- Unique self-closing drawer runners
- High grade 304 stainless steel
- Draws 125mm High
- Shelf & Rails

All stainless steel trolleys are Bulk Ship



Item #	Spec	Dimensions	UOM
1332465	No Drawer	500x500x900mm	Each



Item #	Spec	Dimensions	UOM
1309609	Single Drawer	500x500x900mm	Each



Item #	Spec	Dimensions	UOM
1335464	Single Drawer with Lock	500x500x970mm	Each



Item #	Spec	Dimensions	UOM
1317458	2 Drawers	500x500x900mm	Each



Item #	Spec	Dimensions	UOM
1336096	2 Drawer Side by Side	800x500x900mm	Each



Item #	Spec	Dimensions	UOM
1323024	3 Drawer	500x500x900mm	Each



Item #	Spec	Dimensions	UOM
1332783	4 Drawer	500x500x900mm	Each

OXYGEN TROLLEYS



Oxygen Trolleys

Item #	Suits	UOM
1093194	Size C & D Cylinders Retractable	Each
1228043	Size E Cylinder	Each

STOOLS



Stool Gas Lift Grey w/ Foot Rest

- Adjustable height: 490-630mm
- Cushion diameter: 350mm
- Foam thickness: 60mm
- Base: 560mm

Item #	UOM
1113591	Each



Care Quip X6 Deluxe Mobile Air Bed/Chair w/ Tray & Footrest

- Ideal for dependent clients who are confined to a chair for long periods
- Pressure relieving memory foam, contoured backrest with lumbar and lateral support, angled armrests and adjustable wings and neck pillow ensure maximum comfort and postural support
- Integrated air cell in seat for premium pressure relief
- Tilt in space recline action and high leg rest elevation for reduced fluid retention
- Designed for use with standing machines and transfer hoists
- Large push handle and simple recline mechanism makes movement and operation safe and easy for the carer
- Resilient European four-way stretch vinyl is fluid resistant, fire retardant and vapour permeable with antibacterial and non-allergenic properties

- Seat Width: 500mm
- Back Height: 650mm
- Seat Depth: 500mm
- Overall Width: 750mm
- Arm Height: 180mm
- Max User Weight: 180kg

Item #	Spec	UOM
1451669	Care Quip X6 Deluxe Mobile Air Bed/Chair	Each
1504215	Foot Rest	Each
1504207	Tray	Each



Aspire Mobile Air Chair

- Unique design promotes skin immersion and envelopment, redistributing skin pressure, reducing the risk of pressure injury associated with prolonged sitting or lying
- Swing away backrest wings & armrests. Allows for easy patient side transfers
- Designed with the aged care environment in mind
- Lockable castors with directional control
- Ergonomic attendant grip handles
- Max user weight 180kg

Item #	Spec	UOM
1368516	Aspire Mobile Air Chair	Each
1368532	Meal Tray	Each



Aspire Adjustable Day Chair

- High back orthopaedic chair that offers adjustability, comfort and support
- Contoured backrest
- Depth adjustable
- Latte vinyl colour

Item #	S.W.L	UOM
1393057	160kg	Each




Aspire Waterfall Day Chair

- Segmented and adjustable backrest cushion provides comfort and postural support
- High back orthopaedic chair that offers adjustability, comfort and support
- Padded seats and armrests
- Height adjustable
- Ink vinyl colour

Item #	S.W.L	UOM
1391402	200kg	Each



MATTRESSES



Icon Pressure Care Mattress

- Profile stabilising sides enables ease of care
- Triple layered, premium pressure care core designed to provide optimal pressure care whilst providing support and comfort
- Premium visco grade, pressure sensitive memory upper minimises resident interface pressure by maximising contact surface area
- Water-based evaporative lamination - completely odour free

Optional Extras Include:

- Safety sides - Soft safety sides provide unobtrusive support to help prevent high risk falls residents rolling out of bed
- High frequency sealed and welded cover seams - Provides maximum infection control and facilitates cover cleaning and care



Icon I3100

Item #	Length	Width	Depth	Notes	UOM
1348485	1980mm	880mm	150mm		Each
1372467	1980mm	1050mm	150mm	*suits user up to 280kg King Single	Each
1358103	2030mm	1050mm	150mm	*King Single	Each

Icon IM1000 Bariatric

Item #	Length	Width	Depth	Notes	UOM
1404555	1980mm	1050mm	180mm	*suit users up to 350kg	Each
1444158	1980mm	880mm	150mm	*120kg-280kg	Each
1497359	1960mm	860mm	180mm	*90kg-350kg	Each

All Icon Mattresses are Bulk Ship




BULK SHIPMENT

MATTRESSES WITH PRESSURE REDUCING INSERT



The Forté Assure Evacuation System eliminates the annoyance of traditional mattress evacuation systems through its totally integrated unique design. With all evacuation straps self-contained within inside the mattress cover, there is no safety or damage issue of straps falling loose and becoming entangled with nearby medical equipment or care staff attending to the patient. With all seams High Frequency Sealed and Welded, the Assure Evacuation System

does not compromise on hygiene and infection control. Staff safety and risk of injury is mitigated as the Assure System significantly reduces force required to evacuate a patient across a range of floor surfaces. Unique foam core modifications allow the patient to be encapsulated and further protected against harm.



Mattress with Customised Roho® Insert

- Designed as a high level product for pressure reduction
- Suitable for use with High Risk and Palliative Care
- Available with any Forte static support surface
- Custom made to resident/facility requirements
- 24 month warranty

Item #	Length	Width	Notes	UOM
1396161			Single Mattress	Each
1297708	860mm	51mm	Roho Insert	Each



BULK SHIPMENT

MATTRESSES WITH PRESSURE REDUCING INSERT



Evacuation Pressure Care Mattress

- Significantly reduces friction forces between mattress and varied floor surfaces to make evacuation procedure easier
- External Seams High Frequency Sealed, Welded and Fused for added infection control
- Wipe Clean Pocket with High Frequency Welded Zipper for Infection Control
- Allows foams to encapsulate patient for protection
- Tapered to allow grip to bed base while in normal use
- Easy to Access, do not need to lift mattress to access
- 4 x Integrated Robust Heavy Duty Carry Handles
- Stainless Steel, Non-Corrosive Nyloc Attachment
- Compliant with Infection Control Requirements
- Fully Integrated & Contained Evacuation Straps
- Prevents Ingress of Fluid through Seams
- 100% Wipe-Clean, Impervious to Fluids
- Pockets Integrated within foam core
- Custom Foam Core Modification
- Low Resistance Base Sheet
- Dual Evacuation Pull-Point

Item #	Length	Width	Depth	Notes	UOM
1454129	1980mm	880mm	150mm	High Risk	Each



Zephair Pressure & Palliative Care Mattress System

High Level Acute Pressure Care and Comfort Mattress System.

The ZEPHAIR Mattress range provides true leadership in powered patient pressure care. Tried and trusted in facilities from ICU Wards to Palliative Hospice use, this pressure care system will provide the answer to all acute care needs, whilst providing unparalleled comfort and staff Primary Care compliance.

Key clinical attributes include ventilated immersion upper layers, deep 24 cell alternating integrated inner mattress core, firm base foam layer, quiet efficient air flow generator with easy set weight calibration and dual stretch waterproof breathable cover.

The ZEPHAIR Pressure Care Mattress range includes models suited for standard hospital bed sizing, wide & bariatric requirements. The quality and features of the Zephair range allow confident usage in the Hospital, Aged Care and Spinal environments. Select the Zephair Pressure Care System and achieve unsurpassed reliability and consistency of care.

24 Cell, Alternating Air System

- 125mm depth facilitates immersion
- Ultra-flexible, soft-touch material
- 24 Cell System minimises patient disturbance
- Static head support region

5 year foam core warranty
2 year air cell warranty
1 year pump warranty



Item #	Spec	Length	Width	Depth	UOM
1395815	Single	1980mm	880mm	215mm	Each
1395823	King Single	1980mm	1050mm	215mm	Each

VITAL SIGNS
MONITORS

WHEELCHAIRS

PRESSURE
MATTRESSES

TILT
BEDS

FIND OUT MORE!



australia HEALTHCARE
SUPPLIES

**ALL YOUR FACILITY SET UP NEEDS
FROM ONE GREAT SUPPLIER**

**VISITOR
CHAIRS**

BEDS

**SYRINGE
DRIVERS**

**LINEN &
PILLOWS**

**PILL
CRUSHERS**

**BEDROOM
FURNITURE**

**NURSE CALL
ALARMS**

SILENT KNIGHT



ALTERNATING MATTRESSES



Smik Care™ Alternating Overlay Pressure Mattress System with Pump

- Assists carers to provide maximum pressure relief
- Multi stretch vapour permeable cover
- Equipped with alternating, static and auto firm (nursing) modes
- Auto firm (nursing) mode maximises the internal pressure
- Retainer air cell design prevents over inflation
- Max user weight: 150kg

Item #	Spec	Length	Width	Depth	UOM
1365657	Single	1950mm	950mm	130mm	Each
1374303	King Single	1950mm	1070mm	130mm	Each
1474081	Single	1980mm	880mm	180mm	Each
1474154	King Single	1980mm	1050mm	180mm	Each
1504436	Long Single	2030mm	880mm	185mm	Each
1504444	Long K/Single	2030mm	1050mm	185mm	Each



Forte Airmonte A5000 Alternating Air Mattress Overlay

- For medium to high risk
- Two-way stretch, vapour permeable cover
- 160kg capacity

Item #	Spec	Length	Width	Depth	UOM
1402862	Single	1960mm	860mm	125mm	Each



Forte Airmonte A8000 Full Replacement Air Mattress System

- High quality pressure redistribution surface suitable for residents with medium to high risk of Pressure Injury
- Specially designed cells incorporated with a low pressure high flow pumping system, enables maximum pressure redistribution

Item #	Spec	Length	Width	Depth	UOM
1402854	Single	1980mm	870mm	150mm	Each



CairMax Duo Foam Air Support Surface System

The CairMax Duo is a versatile foam air support surface with the ability to alternate between static or dynamic support, in response to the needs of the patient.

Mattress functions as a static pressure relieving surface without the control unit, however when clinically required, it can step up to provide dynamic support with the control unit connected.

No concerns over correct pressure level settings with AutoCair automatic weight sensing technology, delivering fail-safe patient comfort and support.

Combining durable foam core with alternating air cells embedded with castellated foam to deliver comfort and support for patients up to 230kg. The castellated foam articulates when the mattress is profiled, while static head area foam offers stable support for enhanced sleep.

Item #	Spec	UOM
1481827	King Single	Each



ALTERNATING MATTRESS



ProCair Prime Mattress Replacement System Pumpless

The ProCair Prime can be used as a fully alternating pressure mattress and can also be easily converted to a hybrid foam mattress.

Easy to set up, the ProCair Prime automatically adjusts the pressure to suit the weight of the patient and does not need a separate pump unit to hang from the bed - it is all built-in! This minimises trip hazards or damage from pumps being bumped off the bed.

No more concerns over correct pressure level settings with our AutoCair. This cyclic automatic weight sensing technology delivers fail-safe patient comfort levels.

The mattress can easily be converted to a hybrid mattress by adjusting the position of a built-in foam layer. Hybrid mode provides even greater comfort for the patient. It easily and quickly converts back into dynamic mode to step-up to complete pressure off-loading if required.

Item #	Spec	Length	Width	Depth	UOM
1481800	King Single	2000mm	1050mm	180mm	Each



ProCair Plus Mattress Replacement System Sealed Base

The ProCair Plus is our premium solution for high-risk patients and offers optimal patient therapy and comfort; fully automated operational features; advanced safety features including back-up battery; and advanced error and service indicators.

No more concerns over correct pressure level settings with our AutoCair. This feature's continuous automatic weight sensing technology that delivers fail-safe patient comfort levels.

Increased protection and therapy for vulnerable heels with an independent heel zone made up of narrower cells, that can be disconnected as required. An inbuilt BackUp Battery unit provides the reassurance of up to 5 hours emergency backup power in event of a power failure or accidental power disconnection.

Fowler Boost mode automatically adjusts mattress pressure to provide increased support under the sacral area when the backrest is raised.

Be alerted when your patient leaves the bed with our inbuilt CairAlert feature. This has an optional warning alarm signal to indicate that the patient has got off the mattress. Harnessing the best in technical innovation, our feature-rich products deliver a complete range of alternating pressure mattress solutions for high-risk patients, ensuring optimal patient therapy and comfort.

Item #	Spec	Length	Width	Depth	UOM
1481819	King Single	2000mm	1050mm	200mm	Each


MATTRESS UNDERLAY



Alternating Mattress Underlay

Item #	Spec	Length	Width	Depth	UOM
1271948	Single	1980mm	880mm	50mm	Each
1376942	King Single	1980mm	1050mm	100mm	Each

PRESSURE CUSHIONS



EquaGel® Cushions

- EquaGel is a dry-polymer gel, it has unbelievable tensile and compression strength and incredible durability
- 2 Layer Gel technology- top layer is a thin wall with tight grid for comfort, the bottom layer is thicker for deep posture support
- Contoured pelvic area of the cushion for additional equalizing pressure distribution
- Black cushion cover is fluid resistant and vapour permeable
- Grid design spreads & reduces pressure over a wider area
- Advanced Protector Cushion - 2.5" tall maintaining 1 3/4" gel below pressure points
- Moderate General Cushion - 2" tall maintaining 1 1/4" gel below pressure points

Moderate General Cushion

Item #	Length	Width	UOM
1327712	16"	18"	Each
1369709	20"	16"	Each
1334360	18"	18"	Each

Advanced Protector Cushion

Item #	Length	Width	UOM
1327720	16"	18"	Each
1334379	18"	18"	Each
1418203	18"	16"	Each
1418157	16"	16"	Each

Adj Advanced Protector Cushion

Item #	Length	Width	UOM
1369717	18"	16"	Each
1381601	22"	18"	Each



Roho® Quadro Alternating Cushion

- ISOFLOR® Memory Control offers shape fitting capabilities while the user is seated
- Facilitates blood flow by fitting, matching and tracking the shape of a person
- Quick and easy on-demand adjustment to maximize function
- Cell size: 90mmx90mm


Item #	Spec	Length	Width	Depth	UOM
1286048	HP Cushion	410mm	410mm	90mm	Each
1328018	Cushion cover for Roho Quadro				Each



Ultimate Bliss Premium Cushion

- A combination of high density and memory foam helps provide pressure reduction and comfort for the at risk resident
- Cover features a vapour permeable and two-way stretch waterproof top with a non-slip base, entire cover can be cleaned for easy maintenance


Item #	Spec	Length	Width	Depth	UOM
1331566	Cushion	450mm	450mm	30mm	Each
1272201	Cushion	450mm	450mm	90mm	Each
1272210	Wheelchair Cushion	400mm	400mm	70mm	Each



Forte Sonnet 2 Core Cushion

- High performance medical grade foams to provide sustained patient care and comfort

Item #	Spec	Length	Width	Depth	UOM
1418939	Cushion	460mm	400mm	70mm	Each
1418920	Cushion	460mm	460mm	70mm	Each
1419455	Pressure relieving heel cushion				Each



DermaSaver®

- MicroSpring textile for great protection
- Reduces friction and pressure against the skin
- Many designs available – ensures an ideal product to suit your residents' needs

Forearm Tube

Item #	Size	UOM
1388088	18 - 25cm	Each
1388096	25 - 30cm	Each
1388118	30 - 38cm	Each
1388126	38 - 46cm	Each

Heel Protector with Toe Cover

Item #	Size	UOM
1376608	28 - 30cm	Each
1387952	30 - 35cm	Each
1387960	35 - 39m	Each
1387979	Over 39cm	Each

Heel Protector Stay Put

Item #	Size	UOM
1325531	28 - 30cm	Each
1325523	30 - 37cm	Each
1325515	37 - 39m	Each
1325507	Over 39cm	Each

Full Leg Tube

Item #	Size	UOM
1325752	18 - 23cm	Each
1325744	23 - 30cm	Each
1325736	30 - 38cm	Each
1325728	38 - 48cm	Each
1325701	48 - 58cm	Each

Shin Tube

Item #	Size	UOM
1387901	18 - 23cm	Each
1387928	23 - 30cm	Each
1336452	30 - 38cm	Each
1387936	38 - 48cm	Each
1387944	48 - 58cm	Each

Palm Pillow

Item #	Size	UOM
1388142	Small	Each
1388150	Large	Each

Knuckle Protector

Item #	Size	UOM
1388169	>20cm	Each
1388177	<20cm	Each

Toe Tubes

Item #	Spec	UOM
1388185	3 x S, 5 x M, 2 x L	P10

Pillow Case P.I. Reduction

Item #	Size	UOM
1388134	66 - 51cm	Each

Oxygen/Feeding Tube Cover

Item #	Spec	UOM
1325760	Tube Covers	P10

Finger Separator

Item #	Size	UOM
1325493	One Size Fits All	Each

Ankle Bumper

Item #	Size	UOM
1377159	23 - 28cm	Each



CareProtect Pedi Boot

- A unique device to distribute pressure for residents with a pressure injury risk
- High density foam provides low load-bearing effect
- Open toe allows foot to breathe
- Seams on exterior for less irritation
- Reusable and washable

Item #	Size	Colour	UOM
1319671	27-31cm	Green	Pair
1319698	32-36cm	Blue	Pair
1319701	37-42cm	Yellow	Pair



Gelbo Elbow / Heel Protector

- Mineral oil-based gel formula: softens and moisturises dry cracked skin
- Reduces pressure, friction and shear forces

Item #	Size	UOM
1240324	Extra Large	Pair

BODY PROTECTION



Care Quip Silicone Fiber Heel Protector

- Padded heel cushion, provides comfort and assists in the prevention of pressure build-up in the heel area
- Blue, one size fits most
- Machine-washable

Item #	UOM
1180772	Each
1383353	Each



Protect-A-Limb®

- Prevents skin tears and damages
- Protection against friction to vulnerable limbs

Item #	Size	UOM
1268989	Small	Pair
1268998	Medium	Pair
1269019	Large	Pair



Sheepskin Heel and Elbow Protector

- Hi-temp washable
- Adjustable velcro fitting
- Green, one size fits most

Item #	UOM
1184190	Pair



Comfy Ears® Large

- Cushioned protection for Oxygen Cannula Tubing
- Can be trimmed to required size

Item #	UOM
1226209	Each



Heelmedix Standard Heel Protector w/ Wedge

Item #	UOM
1375814	Each



Pressure Relieving Heel Wedge

Item #	Size	UOM
1225640	43 x 22 x 9cm	Each
1216033	86 x 46 x 9cm	Each
1215878	45 x 18 x 9cm	Each



Pressure Relieving Body Wedge

Item #	Size	Height	UOM
1217729	60 x 30cm	9 - 14.5cm	Each



Molnlycke Z Flo Fluidised Flex Heel Boot

Fluidised heel boot for offloading pressure to help prevent and manage heel pressure ulcers.

Item #	UOM
1415859	Each



Molnlycke Z Flo Fluidised Med Positioner w/ Tab

Easily-customisable fluidised pressure redistribution device to assist in the safe turning, positioning and repositioning of adult patients

Item #	Size	UOM
1415867	41 x 76cm	Each



Molnlycke AD Large Positioner

Molnlycke® Z-Flo™ Fluidized Positioners assist in turning, positioning and repositioning adult and geriatric patients safely and effectively.

Item #	Size	UOM
1484257	63.5 x 91.5cm	Each



HipSaver Nursing Home

Designed for the mobility and comfort of residents plus convenience and peace of mind for caregivers. High Compliance Nursing Home Model has permanently sewn-in protective pads - for higher resident compliance.

Protection that's easy to wear & manage. Designed for the mobility and comfort of residents plus convenience and peace of mind for carers. HipSaver Nursing Home is a durable brief with a roomy fit for easy dressing and continence care. All HipSaver garments and pads are laundries safe for washing and drying up to 95°C.

- Durable, breathable polycotton-spandex brief
- Ample leg openings for easy dressing and toileting
- Soft, thin protective hip pads
- Fits easily over diapers, accommodates incontinence pads
- Can be worn 24 hours a day - comfortable to sleep in

Starter Kit

Item #	Size	UOM
1336436	Medium (F)	Each
1336444	Medium (M)	Each
1395688	Small (M)	Each

High Compliance

Item #	Size	UOM
1333666	Small (F)	Each
1335847	Medium (M)	Each
1406809	Large (M)	Each

Pant Only

Item #	Size	UOM
1325132	Large (F)	Each
1325140	Medium (F)	Each
1325159	Small (F)	Each
1325191	Extra Large (M)	Each
1325205	Large (M)	Each
1325213	Medium (M)	Each



SlimFit Pant Only

HipSaver SlimFit delivers excellent protection. Made from polycotton with form fitting Spandex®, HipSaver SlimFit looks good, feels good and pulls on just like regular underwear. It's lightweight and flexible enough for walking, fold, tennis, skiing and other sports, and comfortable enough for everyday wear.

HipSaver SlimFit is ideal for fashion conscious and active elders who are at risk of hip fracture, but reject uncomfortable and obvious hip protection products.

- Soft, thin (13mm) pads, tapered for protection with style
- Ultra-soft breathable polycotton-spandex for comfort and smooth fit without being tight or binding
- Comfortable enough to wear 24 hours a day - even when sleeping
- Fully machine washable and dryable, even at high temperatures
- Sizes from XS to XXL for a proper fit

Item #	Size	UOM
1325264	Extra Large (F)	Each
1325272	Large (F)	Each
1325280	Medium (F)	Each
1325299	Small (F)	Each
1325345	Large (M)	Each

Hip Protecting Pads for all HipSaver Garments



Item #	UOM
1325477	Pair

SHEETS & PILLOWCASES



PolyCotton Stretch Knit Fitted Sheets

- Designed to encase the mattress and overlay
- Allows the body to conform to the pressure reduction features of the support surface - without hammocking

Item #	Size	UOM
1263721	900x2000mm	Each



Excel Flat Sheets & Pillowcases

- Made to Australian standards
- Will not pull
- Only 2% shrinkage
- 50% polyester / 50% cotton


Item #	Spec	Size	UOM
1274912	Pillowcase	780x500x140mm	B100

BLANKETS



Ultra Plush Microfibre Antished

Item #	Colour	Thickness	UOM
1419544	Charcoal	300gsm	Each
1419552	Mushroom	300gsm	Each



Coral Fleece Large Single

Item #	Colour	Size	UOM
1391410	Walnut	1800x2450mm	Each
1391437	Chamois	1800x2450mm	Each



Celloweave Blanket

- 100% cotton preshrunk white

Item #	Colour	Size	UOM
1337025	Single	1800x2300mm	Each

BATH TOWELS & FACE WASHERS



Pearl Bathing Range

- Designed to withstand commercial laundries and perform for an extended life
- 100% cotton carded pile, soft to the touch and highly absorbent

Item #	Spec	Colour	Size	Thickness	UOM
1308890	Bath Towels	White	1550x710mm	550gsm	C48
1291092	Face Washers	White	330x350mm	56gsm	P50



Other Bath Towels & Face Washers

Item #	Spec	Colour	Size	Thickness	UOM
1290908	Face Washers	White	280x280mm	Small	Each
1338811	Face Washers	White	330x330mm	55gsm	Each
1316443	Face Washers	Navy	330x330mm	50gsm	C100



Vinyl Pillow Protector

Item #	Spec	UOM
1029522	Plastic Zippered	Each



Pillow WipeClean

Item #	Spec	UOM
1029522	Waterproof 700mm x450mm	Each

LINEN PROTECTION



Pinkie Pads

- Waterproof
- Replaces macintoshes and drawsheets
- Large size Pinkie Pads can absorb up to 1L of fluid – with or without tuck-ins
- Small size is perfect to be used as a chair pad or a pillow protector

Item #	Spec	Size	UOM
1161721	Large	900x850mm	Each
1287346	Large	900x850mm	Each



Macintoshes & Drawsheets

- Heavy duty and high quality makes them commercial laundry compatible
- Draw Mac can be used in conjunction with Draw Sheets
- Draw Sheets have a herringbone twill for greater absorbency

Item #	Spec	Size	UOM
1218987	Draw Mac - White	1750x900mm	Each
1289632	Draw Mac - Vintex Full Wrap Tuck-in	1800x900mm	Each
1275956	Draw Mac - Heavy Duty with Wings	1000x900mm	Each
1290363	Draw Sheet - Herringbone	1800x1200mm	Each



Absorbent Pads

Item #	Spec	Size	UOM
1286528	Bed Pad - 2L Capacity with Wings Non-waterproof	1000x1000mm	Each
1104189	Chair Pad - Small 300ml Capacity Tartan	450x500mm	Each
1316877	Chair Pad - Large 1L Capacity Dark Blue	500x600mm	Each
1330969	Chair Pad - Large 1L Capacity Dark Green	500x600mm	Each



Drycare Laminated Linen Protector

- Features non-slip waterproof backing

Item #	Spec	Size	UOM
1290916	1L, Waterproof	850x900mm	Each



Plastic Sheeting

Item #	Spec	Size	UOM
1253752	Plastic waterproof sheet	900x1200mm	P100

LAUNDRY BAGS



Reusable Cloth Laundry Bags

Long-lasting heavy duty polyester twill

Standard

Item #	Colour	Size	UOM
1274752	Blue	750x700mm	Each
1274805	Green	750x700mm	Each
1276407	Purple	750x700mm	Each
1274761	Red	750x700mm	Each
1274770	White	750x700mm	Each
1274799	Yellow	750x700mm	Each

Cytotoxic

Item #	Colour	Size	UOM
1281208	Purple	750x700mm	Each



Impermeable Laundry Bags

- PVC coated polyester
- Waterproof and leakproof for greater infection control
- No need for pre-bagging – bag and contents are washed together on foul cycle
- Effective for approximately 200 washes
- Drying temperature is 75°C

Item #	Colour	Size	UOM
1274500	Red	750x850mm	Each
1274529	Yellow	750x850mm	Each



Limiting Sling

- Designed to prevent overfilling
- Results in better washing practices
- Will last for over 500 washes

Item #	Colour	UOM
1278215	Blue	Each
1277146	Green	Each
1295039	Purple	Each
1277138	Red	Each
1276840	White	Each
1278207	Yellow	Each

Short Laundry Bags

Item #	Colour	Size	UOM
1348965	Blue	750x700mm	Each
1335561	White	750x700mm	Each



Soluble Bag Liners

- Designed to fit into most laundry bags as an impermeable liner
- Manufactured from polyethylene

Soluble Seam Liners

Item #	Colour	Size	UOM
1408054	Purple	720x990mm	C250
1274485	Red	720x990mm	C250
1274494	Yellow	720x990mm	C250

Fully Soluble Liners

Item #	Colour	Size	UOM
1286501	Clear	666x838mm	C200



Mobile Linen Trolleys

- Compact stainless steel tubular frame
- Foot operated lid for ease of use
- Durable acrylic lid

Item #	Spec	UOM
1323237	Single Skip Trolley	Each
1323253	Double Skip Trolley	Each
1323261	Triple Skip Trolley	Each
1434373	Clean Linen Trolley Front Cover	Each
1434403	Clean Linen Trolley Large 3 Shelf	Each



Colour ID Strips for Linen Trolleys



- Allows you to colour code linen trolleys for infection control

Item #	Colour	UOM
1371118	Blue	Each
1371126	Green	Each
1411551	Red	Each
1411578	Yellow	Each



Advise colours required

Item #	Spec	UOM
1374710	Custom Spray Paint for Linen Trolley Lid	Each



ONLINE ORDERING MADE EASY!

SSS Australia wants to make ordering your healthcare supplies quick, easy, and cost effective. Our website is designed with you in mind, offering many tools to help you spend less time ordering, and more time providing care. Here are just some of the features of our industry leading website.

1. Greater functionality to manage how you order including flexible templates, order history and custom catalogues
2. More control over your account, including user and site management, budget and approved product protocols, and more
3. View current orders status and all previous orders from SSS Australia
4. Generate order history and back order reports
5. Robust product search by name, brand, item code, vendor code, and more

SYRINGE PUMPS



BD BodyGuard T Syringe Pump

The BD BodyGuard™ T syringe pump is a small, lightweight ambulatory infusion pump which offers flexibility and reliability. Effortless transition between hospital and home care requirements is the core of the BD BodyGuard™ T syringe pump.

Item #	Spec	UOM
1453157	BD BodyGuard T Syringe Pump	Each
1451316	9V Battery Model 6LR61 for BodyGuard T	B10
1264077	Lockbox for BodyGuard T	Each
1330381	Label Subcutaneous Use 60mm x 50mm Brown	Roll
1309021	Label Line Subcutaneous Brown	Roll
1015349	I.V. Admin Extension Tubing 150cm Luer Lock	Each
1264068	Lockbox Carry Bag for Niki & BodyGuard T	Each

MEDICATION ROUNDS



Silent Knight Mechanical Pill Crusher

- Simple to operate – reduces repetitive injury
- Pouches contain tablet during crushing – no mess
- Portable and sturdy – easy to store on medication trolleys

Item #	Spec	UOM
1240914	Silent Knight Precision Pill Crushing Machine	Each
1240941	Silent Knight Pill Crusher Pouches Deluxe	Pack 1000



Pill Cutter & Counter

- Built in magnifier that helps you align your pill or tablet
- Safety shield completely covers the cutting blade when medication is split

Item #	Spec	UOM
1365088	Pill Cutter Ezy Cut with Safety Shield & Magnifier	Each
1368052	Triangular Pill Counter	Each



Nurses Utility Pouch

Item #	Size	Colour	UOM
1435337	160x200mm	Violet	Each
1435345	160x200mm	Black	Each
1435353	160x200mm	Pink	Each
1435418	160x200mm	Blue	Each



Porcelain Mortar & Pestle

Item #	Size	UOM
1072336	100ml	Each



Do Not Disturb Medication Apron

Item #	UOM
1303996	Each

BLOOD PRESSURE - AUTOMATIC



Omron HEM-7156T Plus Blood Pressure Monitor

Accuracy is one of the most important requirements in blood pressure monitoring. With IntelliWrap™ 360° accuracy, OMRON's HEM-7156 allows you to get accurate readings effortlessly regardless of the angle the cuff is positioned on the upper arm. Get the same accuracy as you would in your doctor's office! This model comes with Bluetooth connectivity.

- Omron Connect App makes it easy to save, track and manage your health data in one place
- Stores 60 measurements in memory with date and time
- Average for last 3 readings within 10 minutes
- Irregular heartbeat detection
- Body movement detection
- Cuff wrapping guide
- Hypertension indicator
- Intellisense technology
- Low battery indication
- Bluetooth connectivity

5 year warranty

Item #	1474758	Each
--------	---------	------



Omron HEM-907 Clinical BP Monitor Kit

- Dual measurement mode
- Environment-friendly, mercury-free
- Easy for operation and training
- Quiet, ideal for the clinical environment

Kit includes:

- Blood Pressure Monitor
- Small, Medium and Large cuffs and bladders
- Cuff tube 1-metre length
- Rechargeable battery pack
- AC Adapter / battery charger

Item #	1203910	Each
--------	---------	------

2 year warranty



Omron HBP-1320 Professional Monitor w/ 2 Cuff/Bldr

- With a rugged design and a practical carrying handle, HBP-1320 is designed for use in professional setting and easily carried in everyday use. The device has been clinically proven to deliver fast and reliable results.
- Designed for routine blood pressure measurements. Enhanced usability for nurses carrying the device from bed to bed and continuous blood pressure measurements.
- Sturdy enough to survive accidental falls, and allow you to carry with ease.

5 year warranty

Item #	1424335	Each
--------	---------	------



Welch Allyn ProBP 2000 Digital BP Device

- Blood pressure readings are captured on inflation in about 20 seconds
- Meets AAMI standards for blood pressure accuracy
- Compatible with Welch Allyn FlexiPort Reusable Blood Pressure Cuffs, from child size 9 up to thigh size 13 (15 - 55 cm range)
- Device includes Adult size 11 FlexiPort Reusable Blood Pressure cuff

Item #	1407821	Each
--------	---------	------

BLOOD PRESSURE - MANUAL



Heine Gamma G5 Sphyg

- Latex-free
- Exclusive pressure-distribution system
- Pressure-proof manometer ensures durability and reliability
- Micro filter protects valve and manometer

Item # 1224901 Each



Heine Gamma G7 Sphyg

- Shockproof. Latex-Free. Upper housing made of aluminium alloy with a shock absorbing ring protecting indicator and glass
- Complete with adult cuff, in zipper pouch
- Large Insufflation Bulb. 50 % larger bulb volume for fast inflation of the cuff

Item # 1224910 Each



Welch Allyn Durashock DS66

- Gear-free DuraShock™ technology
- Certified accuracy to +/- 3mmHg
- Unsurpassed reliability with a unique trigger-release design
- Withstands a 30" drop

Item # 1146251 Each



Sphyg Complete Aneroid

- Measurement Range: 0-300mmHg
- Cuff Material: Nylon * Scale Glow-in-Dark *

Item # 1036545 Each



Spirit Aneroid Sphyg

- Large Dial for easy reading of measurements
- Latex free
- No pin stop calibrated scale

Item # 1178845 Each



Welch Allyn 767 Sphyg

Lifetime calibration warranty reduces cost of ownership

10 year Warranty

Item # 1222700 Each



Welch Allyn Flexiport Cuff Range

Affordable and easy to use, Welch Allyn FlexiPort Reusable Blood Pressure Cuffs offer all the advantages of FlexiPort technology - including simplified connectivity - in a dependable, durable design.

Item #	Size	UOM
1287281	10 (200-260mm)	Each
1286110	11 (250-340mm)	Each
1275689	11 (250-340mm)	Each
1286129	12 (320-430mm)	Each
1274171	Long 12 (320-430mm)	Each
1274162	Thigh 13 (400-550mm)	Each



Welch Allyn EcoCuff

Traditional single-patient-use cuff solutions can be expensive, wasteful, and difficult to implement. EcoCuff design gives you the benefits of switching to single-patient-use, including reduced cross-contamination while controlling costs and the volume of hospital waste.

Item #	Size	UOM
1439367	10 Small Adult (200-280mm)	C100
1439359	11 Adult (270-380mm)	C100
1471457	12 Large Adult (330-450mm)	C100

VITAL SIGNS MONITORING



Welch Allyn Connex Spot Vital Signs Monitor

The Connex® Spot Monitor features an easy-to-use, vivid touchscreen display and provides accurate vital signs measurement including blood pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. The device connects wirelessly to your EMR and sends vitals to the patient chart from the point of care to help improve efficiency. Its upgradeable design and custom configurations include optional wall unit integration, industry-leading thermometry and pulse oximetry choices, connectivity and other mounting/mobility options. In addition, the Accessory Power Management Stand with WhisperDrive™ technology features additional power backup, providing up to 17 hours of on-time, alleviating the need to charge the monitor as frequently.

- Customizable workflow options for documenting patient observations (ie. Is & Os, Pain Scale) and vitals modifiers (ie. BP cuff location)
- Automated Early Warning Scoring to help clinicians identify and respond to signs of patient deterioration (ie. MEWS, NEWS)
- Options for mobile stands, wall mounts and integrated exam systems to solve workflow requirements
- Single sign-on options—including integration with Imprivata—available for quick, secure clinician login
- Secure: FIPS-approved encryption of transmitted data and a range of additional security features
- SureTemp® Plus Oral / Axillary or Braun ThermoScan® PRO 6000 Tympanic Thermometry
- Options for WiFi®, Ethernet and Bluetooth® connectivity to EMR systems
- Choice of leading SpO2 technologies: Masimo®, Nellcor® or Nonin®
- BP averaging to help ensure accurate hypertension diagnosis
- Connectivity to Detecto, Health o meter® and SECA scales
- SureBP® technology can provide BP in only 15 seconds

Item #	Spec	UOM
1347101	Connex Vital Signs Monitor NIBP/Nonin SpO2/Pro6000	Each
1348361	Mobile Basket Stand	Each
1396501	Mobile Work Stand	Each



Welch Allyn Spot Vital Signs 4400

The Welch Allyn Spot Vital Signs 4400 Device gives you a simple, all-in-one, touchscreen solution for measuring patient vital signs. Now you can quickly capture a full set of vitals - including temperature, blood pressure, pulse rate and SpO2 - right from the home screen.

- Options for documenting patient observations (like respiration rate and BMI)
- EMR connectivity helps your care team avoid documentation errors
- SureBP technology can record blood pressure in just 15 seconds
- BP averaging helps ensure accurate hypertension detection
- 7" color touchscreen for a simplified workflow
- SureTemp Plus oral/axillary thermometry
- Nonin® SpO2 technology

Item #	Spec	UOM
1463756	Vital Signs Monitor 4400 w/NIBP/Suretemp/Nonin/Plug 6	Each
1463640	Mobile Stand	Each



OTOSCOPES



Welch Allyn Pocket LED Diagnostic Sets

The Welch Allyn Pocket LED diagnostic set combines function and style whilst providing maximum performance and durability in a lightweight, compact set. It is the first pocket size LED diagnostic set that uses advanced SureColour LED technology, which delivers true tissue colour rendition.

Item #	Spec	UOM
1342738	Set 2 Handles in Soft Case Blue	Each
1342703	Set 2 Handles in Soft Case Black	Each
1412159	Otoscope Only Blue No Case	Each
1342967	Otoscope Only in Soft Case Blue	Each



OTOSCOPIES



Heine Mini3000 Otoscope

Modern pocket otoscope in a unique, compact design. Handy and powerful.

- A tip for a better future. HEINE AllSpec disposable tips – now made of recycled plastic
- Easy instrumentation. Viewing window fixed and can be swivelled to the side
- Pneumatic eardrum testing by means of an integrated insufflation port
- Switches off automatically as soon as it is clipped to the pocket
- Sturdy and non-slip due to a high-quality handle housing
- Always close to hand with practical attachment clip
- High-resolution images, minimal reflection
- Natural colour rendering
- 3x magnification
- Detailed view

Item #	Spec	UOM
1257161	2.5V XHL Otoscope w/ Handle & Tips NO Battery	Each
1443518	LED F.O. Otoscope	Each
1467344	LED F.O. Otoscope w/ Hard Care	Each

THERMOMETERS



Welch Allyn Braun Thermoscan Pro6000

- 60 second pulse timer assists with the measurement of pulse rate, respiration
- Innovative PerfectTemp technology adjusts for variability in probe placement
- ExacTemp technology detects the stability of the probe during measurement
- Accurate, reliable and comfortable temperature readings

Item #	Spec	UOM
1337564	Thermoscan Pro6000 with Small Cradle	Each
1337580	Charging Dock	Each
1337572	Rechargeable Battery Pack	Each
1337599	Ear Probe Covers	B200



Omron Thermometer Ear Probe Covers TH839S

- Suitable for Omron TH839S Ear Thermometer
- Replacement Probe covers
- Pack of 40
- Includes applicator ring

Item #	Spec	UOM
1307754	Ear Probe Covers	B40



Infrared Forehead Thermometer

- Measure Distance is between 1-3cm
- Hold 20 readings in memory
- Colour coded temperature reading

Item #	Spec	UOM
1431803		Each

PULSE OXIMETER



Aeon Pulse Oximeter Finger Adult LCD Multi Display

Item # 1319124 Each



HeartSure Fingertip Pulse Oximeter

Item # 1339931 Each



Biolight Pulse Oximeter Finger M70C

Item # 1465805 Each

BLADDER SCANNERS



BioCon-900D Bladder Scanner w/ Docking Station

- The BioCon-900 is the latest release from Mcube, offering a wireless, non-invasive 3D ultrasound system that measures the volume of urine in a patient's bladder
- The BioCon-900 uses ultrasound signals to determine the bladder's outline and total volume of urine. The revolutionary real-time Pre-Scan function locates the bladder prior to scanning and improves the accuracy of results
- Immediate printout of results via the internal printer or use the included CubeScan PC software to upload saved data to a computer for later review

12 Month Warranty

Ultrasound Probe:

- Dual frequency (2MHz, 3.4 MHz),
- B-mode scan image,
- Scan angle: 120 deg,
- Penetration depth (normal patient) 23cm

Item # 1406426 Each



Bladder Scanner Touch Screen Printer Kaixin BVT01

The portable bladder scanner BVT01 is a non-invasive ultrasound scanner used to measure bladder volume for assessing urinary retention and urinary incontinence. It aims to reduce the conduction time urinary frequency and reduce the risk of urinary tract infection and for treatment of urinary system diseases.

Non-invasive ultrasound scanner used to measure bladder volume for assessing urinary retention and urinary incontinence. Features smart and compact design with touchscreen, simple operation and advanced calculating method ensures the correct result for every anomaly shape.

Warranty:

- 24 Month main unit
- 12 Months Probes

Item # 1428101 Each

STETHOSCOPES



Littmann Classic III Stethoscopes

- 3M™ Littmann® Classic III™ Monitoring Stethoscope
- Monitor and assess a wide range of patients Detect normal and abnormal sounds and rhythms
 - Useful in non-critical care environments such as a medical office, general ward, OB/GYN, ambulatory clinic or urgent care
 - Tunable, dual-sided stainless steel chestpiece with open or closed bell
 - Non-chill bell sleeve for greater patient comfort
 - Soft-sealing ear-tips provide optimal comfort and excellent sound occlusion
 - Made in the USA without natural rubber latex

5-year warranty

Item #	Colour	UOM
1333518	Lavender	Each
1334050	Black Chestpiece	Each
1337149	Black	Each
1337785	Navy Blue	Each
1337807	Burgundy	Each
1337858	Ceil Blue	Each
1376470	Turquoise	Each



Spirit Deluxe Dual Head Stethoscopes

- Features chest piece made of aluminium material with metal retaining ring, simplex; binaural and 1 piece molded PVC tubing

Item #	Colour	UOM
1203385	Neon Orange	Each
1153069	Pink	Each
1153078	Aqua Blue	Each
1153031	Red	Each
1153040	Burgundy	Each



Spirit Nurses Complete Stethoscope

- Silicone eartips with new snap-on connections for the ultimate in wearing comfort and acoustic seal
- Scope ID Tag included
- Weights 165g
- Overall length 780mm
- Latex-free

Item #	Colour	UOM
1203198	Orange	Each



Liberty Single Head Stethoscopes

A well made standard nurses single head stethoscope, ideal for taking blood pressure and basic auscultation.

- Flexible 560mm PVC tubing in 3 fashionable colours for departmental coding or personal identification
- Light weight aluminium single sided chest piece
- Ideally suited for blood pressure measurement
- Adjustable chrome plated binaurals
- Supplied in clear PVC display case
- Latex-free

Item #	Spec	UOM
1043049	Royal Blue	Each
1043076	Red	Each
1313932	Black	Each



Liberty Dual Head Stethoscopes

Featuring a two-sided chest piece, only a truly extraordinary precision instrument combines outstanding performance with distinctive styling. Acoustic response is remarkable from its quality engineered chest piece. A very economical and well designed unit.

- Light weight aluminium reversible head chest piece
- Adjustable chrome plated internal binaural
- Flexible 560mm PVC tubing in 4 fashionable colours for departmental coding or personal identification
- Rubber non chill ring for patient comfort
- Stethoscope tubing length 762mm
- Packaged in a hanging display case
- Name tag include
- Latex-free

Item #	Spec	UOM
1314009	Complete Black	Each
1316540	Navy	Each
1316486	Teal	Each
1316508	Purple	Each



Healthweigh Wheelchair Scale

The Healthweigh® Wheelchair Scale allows practices to measure patients confined to a wheelchair in a safe and comfortable manner. Patients can be placed onto the scale via the scale's two wide accessibility ramp. The Healthweigh® Wheelchair Scale also includes movement compensation technology, which eliminates involuntary movement during the weighing process to ensure accurate results. The scale also features a manual/automatic tare option, which allows the weight of wheelchair to be reduced in calculations.

Item #	Capacity	Size	Graduation	UOM
1392719	360kg	1250x940x820mm	0.1kg	Each



Weighing Chair

- Full digital calibration (FDC) ensuring accurate and simple calibration
- Safety wheels with brakes, standard foot rest
- Auto power off to preserve battery life
- Comfortable white plastic seat
- Dry battery or AC power

Item #	Capacity	UOM
1049786	150kg	Each



Seca Weighing Chair 952

Item #	Capacity	UOM
1330691	200kg	Each



Weighing Chair Scale

Item #	Capacity	UOM
1360914	300kg	Each



EXAM LIGHTS



Welch Allyn GS 300 General Exam Light

- Perfect combination of intense light output and cool operation
- Affordable price and rugged performance
- 20,000 Lux at typical working distance
- Ergonomic power connection
- Single White LED
- Touchless on/off

Item #	Spec	UOM
1280554	GS300 with Mobile Stand	Each
1280562	GS300 with Table/Wall Mount	Each



Heine EL3 LED Examination Light

HEINE has set a new standard which maintains that only the best is good enough, from the selection of materials to processing, from light intensity and thermal management to a colour rendering index as high as possible.

Shadow-free illumination. The compact illumination head (Ø approx. 60 mm) allows convenient, precise positioning, particularly in difficult application situations. Optimal positioning eliminates unwanted shadows in the illuminated field.

Item #	Spec	UOM
1330829	w/ Wall Mount	Each
1337416	w/ Mobile Base	Each
1468898	Wheeled Stand Metal Base	Each



POINT OF CARE TESTING



Clinitek Status+ Urine Chemistry Analyzer

- Simple, point-of-care urinalysis with enhanced clinical information.
- Easy-to-use, automated analyzer designed to deliver quick, accurate, reliable results
- Identifies type of strip automatically
- Checks integrity of each test strip
- Comprehensive urine test menu
- Printed record

Item #	Spec	UOM
1296175	Clinitek Status+ Analyser	Each
1163181	Thermal Paper Rolls	B5
1312502	DCA/Clinitek Printer Labels	B5
1005899	Multistix-10SG Strips 100	Bt100

ANAEMIA SCREENING



HemoCue HbA1c 501 Analyzer

HemoCue® HbA1c 501 System is a fast, easy and reliable point-of-care system for hemoglobin A1c tests. It is a fully automated HbA1c POC test that provides reliable results straight away for efficient diabetes care.

- Stop sending samples to labs and waiting for results
- An automated system for HbA1c testing
- Treat and follow up your patients fast and efficiently
- Advantages of using our Point-Of-Care system
- All-In-One kit to provide your patients with direct feedback
- Onscreen instructions simplify every step of the process
- Room temperature storage gives you less to think about
- Certified according to IFCC and NGSP to ensure quality
- Get the right answers at the right time every time

Item #	Spec	UOM
1453149		Each

GLUCOMETER



Accu-Chek Guide Me Meter MO

- Ideal for medical centres
- Fast, easy-fill strip ' results in less than 4 seconds.
- Automatically turns on when you insert a test strip
- Widest dosing area among leading brands; Apply a small drop of blood anywhere along the yellow edge
- No touch strip removal; the strip ejector button allows for fast and hygienic removal of test strips
- Tiny sample size <0.6 'L with under-dosing detection
- Spill-resistant SmartPack vial lets you take one and spill none
- Simple to see, day or night'strip port light for testing in daytime, nighttime or anytime
- Strips not included with Accu Check Guide (Sold Separately)
- Includes new child resist lock

Item #	Spec	UOM
1467816	Accu-Chek Guide Me Meter MO	Each
1369512	Test Strips x100	B100
1377302	Control Solution Level 1 & 2	Each

STAND ALONE ECG



Welch Allyn® CP150 Interpretive 12 Lead Resting ECG Touch Screen

Stand Alone ECG with Interpretation.

- Battery or mains operated for ultimate flexibility
- Advanced filters for optimal ECG trace quality
- 7" touch-screen display and keyboard
- Three lead ECG preview


Item #	Spec	UOM
1314270	W.A. CP150 ECG	Each
1316516	Hospital Cart w/Out Cable Arm/Shelf	Each
1316524	Cable Arm & Shelf	Each



W.A Cardio Suite ECG AM12 Wired Acquisition Module

Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio AM12 wired acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10)

Item #	1438700	Each
--------	---------	------



Cardioline ECG 200S A4 w/ Glasgow Interpretation and PDF 3

Cardioline ECG 200S Electrocardiograph is an advanced 12-lead ECG, which can display, acquire, printout (on A4 paper 12, 6 or 3 channels) and store ECG exams for adults, children and infants.

ECG 200S extracts the critical ECG parameters and offers interpretation with the Glasgow program available as standard. This interpretation algorithm takes into account specific criteria such as patient sex, age, gender and race, and ensures a high level of accuracy. This interpretation must be reviewed and validated by a physician who will complete the medical report before storage and exportation.

Item #	1446010	Each
--------	---------	------



Red Dot Resting ECG Electrodes

Tab-style electrode has pressure-sensitive conductive adhesive that is easy to apply and remove, and provides rapid and intimate skin contact to ensure high quality trace pickup.

Item #	1009940	P100
--------	---------	------

OXYGEN CONCENTRATOR



Compact 525 5 Litre Oxygen Concentrator

Delivering up to 5 litres per minute with high oxygen concentration across all flow rates and adopting the latest developments in the reduction of sound and power consumption, the latest Compact 525 increases patient comfort and lowers operating costs.

- Intelligent power management system utilises turn-down® technology providing less power consumption below 2 l/min flow rates
- Weighing only 16.3 kg the Compact 525 is the smallest and lightest 5 litre concentrator Visual and audible alarms for low Oxygen levels, power failure, pressure drop and service required from DeVilbiss Healthcare
- Oxygen Sensing Device for continuously monitoring Oxygen levels
- New maintenance free innovative rotary valve
- Recessed humidifier nook prevents damage
- Circuit breaker with reset function

Item #	1442759	Each
--------	---------	------



SUCTION PUMPS



Liberty Suction Pump 7E-A 240V 18L/min

The Liberty 7E-A Portable Suction Unit, is a new design, oil-free lubrication pump and based on extensive experience and development trends from domestic and overseas products.

The Liberty 7E-A Suction Pump is specially designed for extracting thick liquid such as blood, phlegm, etc. Not only can this unit be used as an aspirator in the operating room or emergency room, but also as a health care product for the home.

- Portable, compact small size:
- Square negative pressure dial
- Light-weight, easy to carry
- Oil-free lubrication pump
- High negative pressure
- Built-in bottle cavity
- Volume: 18 Litres
- Large Flux
- Low noise
- Power: 240V

Item #	Spec	UOM
1272954	Liberty Suction Pump	Each
1283499	1L Canister	Each
1297201	Disposable 850cc Canister	Each
1369172	Filter SPBF	Each
1408186	Tubing Set (Canister to Filter)	Each
1418653	Canister Lid Only	Each



Askir 30 Surgical Aspirator/Askir 36BR

- Easily transportable
- Highly heat resistant and electrically insulated plastic
- Overflow valve and vacuum indicator

Item #	Spec	UOM
1272963	Askir Suction Pump 30L 230V (1L Jar)	Each
1323822	Askir Suction Pump 36BR with Built in Battery Pack	Each
1309595	1L Jar (Autoclavable)	Each



Disposable Items for Askir 30/36BR

Item #	Spec	UOM
1439952	Askir Round Spacer	Each
1323814	Flovac Askir Reusable Container 1L	Each
1313150	Flovac Disposable Liners	Each
1439944	Flovac Conical Connector (Yellow)	Each

Replacement Items For Askir 30/36BR

Item #	Spec	UOM
1276354	Antibacterial Filter (also fits Liberty)	Each
1428411	Askir Suction Tubing Set 3 Tubes 6 x 10 x 140cm (1), 6 x 10 x 13cm (2)	Set

NEBULISER THERAPY



AiroMax Nebuliser Therapy Unit

The AiroMax Nebuliser is equipped with a powerful and quiet piston compressor that does not need lubrication, providing a reliable low maintenance nebulising system. It is designed to nebulise medication continuously at home or in clinical settings for children and adults for the inhalation of all prescribed bronchial medication's. The high flow design is manufactured to maintain 12-16L/minute flow rate and a nebulisation rate of 0.50mL/minute. All the AiroMax accessories are phalate-free. The nebulising efficiency of the AiroMax Nebuliser makes it an outstanding choice in administering nebuliser therapy.

- Compressor max pressure 2.5 bar, Power 240V, 50Hz, 1.0 Amp
- Masks, bowl, tubing, filters and carry case
- Allows patient-specific individualisation
- Nebulisation rate 0.55ml/minute
- Quiet piston operation
- Five year warranty
- Noise level 52 dbA

Item #	Spec	UOM
1152348	AiroMax Nebuliser	Each

DEFIBRILLATOR



Defibrillator Lifepak CR2 Essential Semi-Automatic

Sudden cardiac arrest (SCA) can happen to anyone—anywhere. Immediate treatment is vital. A victim's chance of survival dramatically decreases for every minute without treatment. That's why public access to defibrillators is so important. They put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best. Layered design with easy-to-follow, bold graphics. Both trained and untrained AED users clearly know how to begin.

- QUIK-STEP™ Electrode Pads: Peel directly off the base for faster side-by-side placement
- Metronome and CPR Coaching: Sets an effective pace and audibly guides users
- Child Mode: Toggle to Child Mode for reduced energy and CPR guidance appropriate for children
- ClearVoice™ Technology: Detects background noise and adjusts the volume for clear instruction
- Highest Available Energy: Up to 360J for more effective shocks as needed
- LIFEPAK TOUGH™: IP55 rating for challenging environments
- 8-year warranty

Item #	Spec	UOM
1415743	Lifepak CR2 Defibrillator	Each
1419005	Lithium Replacement Battery	Each
1419013	Wall Mount Bracket - White	Each
1418998	Quik-Step Essential Replacement Electrodes	Each

Heartsine Defib Samaritan AED



The HeartSine Samaritan AED 500P is a fully automatic defibrillator, available with or without Wi-Fi connectivity, offering excellent value and robust protection against dust and water, all within one of the smallest and lightest designs on the market.

- User-friendly interface and guidance
- Automatic shock delivery if needed
- Remote Wi-Fi monitoring of device status
- Highly portable and lightweight
- Clinically proven technology
- Exceptional dust and water protection
- Cost-effective solution

Item #	Spec	UOM
1289675	Adult Pad Pak 500P CPR Advisor	Each
1500805	Adult PadPak Comp 360P	Each
1289160	Adult PadPak Comp 350P	Each

ssaustralia HEALTHCARE
SUPPLIES

Does More For you

sales@sssaustralia.com.au

1800 777 518

www.sssaustralia.com.au